



Laboratory Equipment Pty Ltd

INSTRUCTION MANUAL FOR TISSUE FLOAT BATHS

Laboratory Equipment Pty Ltd
"Proudly Australian Owned and Operated."
26 Farr Street, Marrickville NSW 2204
Phone +61 02 95602811 Fax +61 02 95606131
www.labec.com.au

Contents

1. Product Introduction.....	1
2. Safety Instructions.....	1
2.1 Must do.....	1
2.2 Don' t do.....	1
3. Working Environment.....	1
4. Product Overview.....	2
5. Technical Parameters.....	2
6. Feature.....	2
7. Control Panel.....	3
7.1 Display the main page.....	3
7.2 Various parameter settings.....	3
8. Cleaning and maintenance.....	4

1. Product Introduction

This product is a split roaster developed by our company using innovative technology. With a personalized operation panel design, advanced structure and novel shape, it is an ideal instrument for modern pathological sectioning.

2. Safety Instructions

2.1 Must do

- Insert the plug into the special socket of the device. Before turning on the power, make sure that the operating voltage of the device is the same as the power supply voltage.
- The product can only use a power socket with a ground wire. Make sure that the plug is firmly inserted into the power socket.
- If the power cord is damaged, in order to avoid danger, it must be replaced by the manufacturer or its maintenance department or similar professional personnel.
- Place the device on a stable, level surface. Do not put it in high temperature and humidity environment.


2.2 Don' t do

- Do not plug or unplug the power plug with wet hands or press any button on the control panel.
- Do not use the device in an environment with electromagnetic interference.
- Do not modify, disassemble or repair this equipment without permission from the manufacturer.
- Do not wipe the instrument with volatile solvents such as acetone and xylene.

3. Working Environment

- a) Suitable ambient temperature: 5 °C ~ 30 °C
- b) Relative humidity: ≤80%

4. Product Overview



5. Technical Parameters


- a) Dimensions: length * width * height=327mm*340mm*166mm
- b) Voltage And Frequency: AC220V±10%, 50Hz±1Hz
- c) Power: 300VA
- d) Temperature setting range: 0°C ~ 70°C, the minimum graduation value is 1°C
- e) Bath Bowl: 210*168*55(mm), 2.5L

6. Feature

- a) High-definition color touch screen, information is clear at a glance, easy to operate.
- b) The computer accurately controls the temperature, the temperature rises quickly and smoothly, and the error is small.
- c) The timing can be set.
- d) The LED light display at the bottom of the booth allows you to see the tissue slices floating on the water more clearly.

7. Control Panel

7.1 Display the main page



1. Current time display: Click to enter the current time setting page
2. Baking slice temperature setting and related display
3. Timing reminder time setting button: click to enter the timing duration setting
4. Timing start button: click to start timing reminder, click the small horn icon to turn off timing reminder function after alarm
5. Heating status display
6. Heating start and stop button
7. The key to turn on the lighting: click to turn on/off the bath lighting;
8. Function setting key: Click to enter the function setting page

Current date and time settings :

1. Click the date and time area on the main page, and the time setting window will pop up.
2. After entering the current time, click the "Finish" button to complete the time setting.

Timed reminder settings :


1. Click "Time Reminder" on the main page, and the value window will pop up.
2. After entering a reasonable value, click "Confirm" in the value window to complete the duration setting.

Temperature setting :

1. Click "target temperature" in the grill area on the main page, and a numerical window will pop up.
2. After inputting a reasonable value, click "Confirm" in the value window to complete the

setting of the temperature.

7.2 Various parameter settings



Timing boot settings:

1. Click the "Settings" button on the homepage to enter the settings page.
2. Click "Timing boot-Hours" to set the boot hours value.
3. Click "Timing boot-Min" to set the start-up minutes.

Timer shutdown setting :

1. Click the "Settings" button on the homepage to enter the settings page ;
2. Click "Timing off-Hour" to set the shutdown hour value ;
3. Click "Timing off-Min" to set the shutdown minutes ;

Equipment debugging (password) :

1. Click the blank space in the upper right corner of the setting page, and the factory personnel will debug and use the heating accuracy.
2. No special reason is not open to users.

8. Cleaning and maintenance

Note: When cleaning the machine, a mild commercial detergent or soap can be used to clean the dirt.

Do not wipe the instrument with volatile solvents such as acetone and xylene. When cleaning, do not put cleaning liquid into the equipment internal. After cleaning, wipe the surface of the machine with a soft cloth.