

Immersion Circulator

(Ambient +5°C to 100°C)

LABEC Immersion Circulator

is available with either digital set and display of actual operating temperatures or analogue settings. The PID controller provides very precise temperature control and is able to adapt to many different containers. Equipped with additional safety cut-off.

OPTIONS:

- Pump set for external temperature applications.
- S - Optional Safety Cut Off.



SPECIFICATIONS

IMMERSION CIRCULATOR/HEATER

Model	M-ICB-12	M-ICD-12
Temp. Range	up to 100°C	
Temp. Stability	+/-0.3°C to 1°C	+/-0.1°C
Controller	Basic Analog	PID Digital
Readout	Thermometer	Digital °C/°F
Pumping Flow	12 Litre/min	
Immersion	Minimum 6cm, Maximum 155mm	
Heater (W)	1000 Watt	
Over temp. Safety	Adjustable	
Dimensions (mm) W x D x H	200 x 240 x 90	
Weight (kg)	2.3	2.3

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
 email: sales@labec.com.au
 Ph: 02 9560 2811 • Fax: 02 9560 6131
 www.labec.com.au