

Steam Steriliser (+100°C)

LABEC Steam Sterilisers are floor standing and are manufactured in Australia entirely from 304 grade stainless steel.

Fast heat up times create rapid steam production and the steriliser is equipped with a safe opening lid with extended lift arm. The large gabled lid with vent removes the steam from the atmosphere.

The constant level water device prevents water from running low and automatically refills the steriliser. A jacketed double wall construction leaves a cool to touch exterior.



Model SS1

Features include:

- Stainless steel vessel.
- Simmerstat control.
- Gabled lid with vent.
- Constant water level from mains supply.
- Safety cut-out device.
- Rapid steam production.



Model SS2

SPECIFICATIONS

Model	SS1	SS2
Int. Dim. (mm) H x W x D	500 x 400 x 400	150 x 450 x 250
Ext. Dim. (mm) H x W x D	1000 x 500 x 600	300 x 550 x 350
Temp. Min	Ambient +5°C	
Temp. Max	100°C	
Capacity (L)	80	16
Power (W)	3600	2400
Weight (kg)	30	22



Custom Model



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
 email: sales@labec.com.au
 Ph: 02 9560 2811 • Fax: 02 9560 6131
 www.labec.com.au

Serving the Scientific, Medical and Research industries since 1945

