

# Breast Milk Refrigerator (+4°C)

**LABEC** offer 3 models that are a purpose designed breast milk refrigerator. This breast milk refrigerator is designed for maternity wards and hospitals for the safe refrigeration of breast milk.

Model BMR-400  
Glass Door

## Features and benefits:

- Glass door or Solid Door.
- Door lock.
- 3 or 5 shelves.
- Castors (400/600 litre models).
- Australian made.
- Breast milk temperature maintained below 4°C.
- Built in alarms.
- Digital temperature indicator.
- Interior lighting.

## Optional:

- Data logger and software.
- External alarm ready BMS.
- Temperature chart recorder with battery back up.



Wireless Data Logger  
(optional)



Model BMR-130  
Glass Door

## SPECIFICATIONS

Model	BMR-130	BMR-400	BMR-600
Temperature °C	Breast milk temperature maintained below 4°C		
Dimensions (W x D x H) mm	540 x 540 x 820	600 x 600 x 1855 + castors	775 x 775 x 1900 + castors
Capacity (L)	130 Litres	360 Litres	600 Litres
Door	Glass or Solid Door		
Shelves	3 shelves	5 shelves	

Please specify when ordering after model number if you require glass or solid door eg. BMR-400/GD

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd  
email: sales@labec.com.au  
Ph: 02 9560 2811 • Fax: 02 9560 6131  
www.labec.com.au