

# Digestion Block Heater (+400°C)

**LABEC** offers a **Digestion Block Heater** for heavy acid environments with stainless steel exterior and controls at 3 meters length. The digital PID temperature controller controls with  $\pm 1^\circ\text{C}$ . The controller also features an alarm cut-off which is adjustable and must be manually reset if a fault occurs. These units are custom built to size.



## SPECIFICATIONS

Model LDBH-1223

Model	LDBH-1223
Number of holes	12
Depth of holes (mm)	90
Diameter of holes (mm)	23
Temperature Control	Digital PID
Temperature Max	400°C
Power	240V / 15AMPS / 50Hz

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd  
email: sales@labec.com.au  
Ph: 02 9560 2811 • Fax: 02 9560 6131  
www.labec.com.au