

Nail Salon

Ultrasonic Cleaner • Autoclave • Water Distiller

Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC Nail Salon provides your salon with a complete solution and answer to your sanitation and sterilisation requirements.



1 THE FIRST STEP AGAINST INFECTION AND CROSS CONTAMINATION

The **LABEC USC-5L Ultrasonic Cleaner** is designed for dental practice environments. This unit incorporating Lateral Oscillation Technology that sanitises your needles and other instruments ready for autoclaving. The prepared equipment is polished to a fine finish free from stain, dust, oil, grease and other substances.



Lid and Baskets included

Model	USC-5L
Int. Dim. (mm) L x W x H	240 x 140 x 150
Capacity (Litres)	5
Transducer Power (Watts)	200
Operating Frequency (kHz)	40
Material (interior & exterior)	Stainless steel
Heating	75 Watts / 80°C
Digital Timer (Min)	1 - 99

LABEC
Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
email: sales@labec.com.au
Ph: 02 9560 2811 • Fax: 02 9560 6131
www.labec.com.au

Nail Salon

Ultrasonic Cleaner • Autoclave • Water Distiller



2 THE SECOND STEP AGAINST INFECTION AND CROSS CONTAMINATION

The **LABEC AA-17S Bench Top Autoclave** is a purpose built self contained sterilisation unit. Able to sterilise wrapped and unwrapped solid instruments. Australian TGA approved. Conforms to EN10360/2/4 CE, AS4815. With sterilisation temps 105°C and 134°C. Single pre-vacuum preceding sterilisation. Bowie & Dick test confirmation of steam penetration. Stand alone USB chart recorder and USB memory stick capable to record all batch parameters for QA documentation.

Model	AA-17S
Int. Dim. (mm) Dia x Length	250 x 355
Ext. Dim. (mm) L x W x D	590 x 465 x 450
Temp. range	105° to 134°C
Capacity (L)	17
Power (W)	1800
Power Source	10 AMP/ 240 Volt
Controller	LCD Display and 7 step programmable controller.
Printer	Included
Safety Features	Automatic door lock under pressure, safety valve, temperature control and pressure protection system together with a temperature protector mounted on the chamber wall.
Weight (kg)	48

3 THE FINAL STEP AGAINST INFECTION AND CROSS CONTAMINATION

The **LABEC WS1.0 Junior Water Distiller** provides you with high grade distilled water. Use this water for manual cleaning, rinse bath, ultrasonic bath and autoclave water. By not providing pure demineralised water throughout the sanitation and sterilisation process there is a carry over of the salts/minerals and other solids that come through our conventional town water supply. Have you looked inside your kitchen kettle lately? Don't let the brown stain of shame into your system.

Model	WS1.0
Distiller Water	1.0 L/Hour
Storage tank (L)	4
Outside Size (mm) W x H x D	320 x 220 x 385
Weight	4.5 kgs
Voltage / Watts	240V / 750W

DISTILLATE QUALITY

pH: 5.4 to 6.0

Conductivity: 10 to 20 µS/cm