

# Ultra Low Temperature Freezer

## Compact Models (-10°C to -40°C)

**LABEC** offer two compact models in our **Ultra Low Temperature Freezer** suitable for under the bench or small spaces. Up to seven draw capacity, Ideal for the storage of small items where limited space is available.

### Control System

**Microprocessor temperature controller,**

**Temp. range:** from -10°C to -40°C. LCD display clearly shows the temperature and functions of the controller.

**Perfect audible/visual alarm:** High or low temperature alarm, system failure alarm. **Power supply:** 240V/50Hz 1 phase.



### Structure Design

- Upright type.
- Moulded plastic material interior.
  - Outer body is painted steel.
  - Door safety lock.
  - Door handle.
- Castors fitted for easy mobility.
- Interior drawers made of plastic.

**Optional:** Chart recorder.



### Refrigeration System

Danfoss compressor and EBM fan motor from Germany. Refrigerant as R507, CFC free.

**Certificate:** CE, ISO9001, ISO14001, ISO13485

## SPECIFICATIONS

Model	Capacity (L)	Drawers	External size (WxDxH) (mm)	Internal size (WxDxH) (mm)	Input Power (W)	Weight (kg)
DW-FL90 (CE MARK)	90	4	623 x 544 x 842	438 x 388 x 539	280	45
DW-FL270 (CE MARK)	270	7	700 x 600 x 1750	440 x 400 x 1490	370	65

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



**Laboratory Equipment Pty Ltd**

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au