

Ultra Low Temperature Freezer

(-15~-40°C / -25~-60°C / -40~-86°C)

Serving the Scientific, Medical and Research industries since 1945

LABEC



H-DW-40L58



H-DW-40L408



H-DW-40L938

Temperature Control Function

- Micro-computer control, digital temperature display, accurate to 0.1°C.
- Temperature is adjustable, temperature difference of start-up and shut-down can be set.
- High and low temperature alarm function, alarm temperature can be set as needed.

Safety System

- Perfect alarm system with buzzer sound and lamp flashing alarm functions.
- High and low temperature alarm and sensor fault alarm functions are available.
- Multiple-protection function - password protection, frequent start-up protection, refrigerating system operation protection.

Refrigeration System

- Ultra-thick foam layer and dual-sealing door for effective temperature protection.
- Optimized auto-cascade refrigerating system with high quality compressor for excellent refrigerating performance.
- Evaporator design to ensure the most effective flow of refrigerating fluid and uniformity of temperature in the cabinet.
- Unique fluorine-free mixed refrigerant, 100% CFC free.
- High quality noiseless fan for stable performance and lower noise.

Design Features

- Food-grade stainless steel liner prevents rust.
- Heavy duty casters for convenient movement.
- Safety door design to prevent accidental opening.
- LED temperature display.



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au

SPECIFICATIONS

Serving the Scientific, Medical and Research industries since 1945

Model	External Dimensions (WDH) mm	Internal Dimensions (WDH) mm	Temperature (°C)	Power (W)	Net weight (Kg)	Capacity (L)	Shelves
DW-40L028	470*490*1120	270*265*400	-15~-40	350	64	28	1
DW-60L028			-25~-60				
DW-86L028			-40~-86				
DW-40L058	540*620*1140	364*405*409	-15~-40	450	73	58	1
DW-60L058			-25~-60				
DW-86L058			-40~-86				
DW-40L108	765*730*1260	460*485*510	-15~-40	550	113	108	1
DW-60L108			-25~-60				
DW-86L108			-40~-86				
DW-40L158	700*700*1580	460*455*800	-15~-40	630	130	158	2
DW-60L158			-25~-60				
DW-86L158			-40~-86				
DW-40L208	875*730*1480	570*800*485	-15~-40	650	138	208	2
DW-60L208			-25~-60				
DW-86L208			-40~-86				
DW-40L308	875*730*1820	570*485*1140	-15~-40	690	149	308	2
DW-60L308			-25~-60				
DW-86L308			-40~-86				
DW-40L408	850*930*1950	550*630*1140	-15~-40	750	230	408	3
DW-60L408			-25~-60				
DW-86L408			-40~-86				
DW-40L508	875*730*1900	700*655*1150	-15~-40	1400	225	508	3
DW-60L508			-25~-60				
DW-86L508			-40~-86				
DW-40L608	940*950*1920	700*680*1255	-15~-40	1600	280	608	3
DW-60L608			-25~-60				
DW-86L608			-40~-86				
DW-40L708	1155*980*1860	850*735*1170	-15~-40	1630	290	708	3
DW-60L708			-25~-60				
DW-86L708			-40~-86				
DW-40L808	1155*980*2025	850*735*1335	-15~-40	1650	300	808	3
DW-60L808			-25~-60				
DW-86L808			-40~-86				
DW-40L938	1240*1000*2000	1000*700*1335	-15~-40	1700	320	938	3
DW-60L938			-25~-60				
DW-86L938			-40~-86				

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
 email: sales@labec.com.au
 Ph: 02 9560 2811 • Fax: 02 9560 6131
 www.labec.com.au