

Ultra Low Temperature Freezer

Small Capacity Chest (-40°C to -86°C)

LABEC offer a range of **Small Capacity Ultra Low Temperature Chest Freezers** ideal for the storage of small items and where limited space is available.

Features include:

- Digital Control system with over and under temperature alarms.
- Anti corrosion stainless steel liner.
- Stable temperature with double-sealed lid gasket.
- CFC-free refrigerant system.
- Washable condenser filter.
- Powder coat steel finish.



SPECIFICATIONS

Model	H-DW-86W58	H-DW-86W118	O-DFC84
Int. Dim. (mm) W x D x H	410 x 410 x 355	540 x 400 x 560	450 x 400 x 470
Ext. Dim. (mm) W x D x H	670 x 760 x 1025	1270 x 780 x 965	700 x 650 x 1080
Temp. Min	-86°C		
Temp. Max	-40°C		
Capacity (L)	58	118	84
Refrigerant	CFC Free		
Controller	Microprocessor Control System		
Power	240V / 50Hz		
Weight (kg)	70	110	150

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
 email: sales@labec.com.au
 Ph: 02 9560 2811 • Fax: 02 9560 6131
 www.labec.com.au