

Desiccator Cabinets

Stainless Steel

LABEC Desiccator Cabinets are constructed of quality stainless steel with a full acrylic door.

All joints are sealed and the adjustable door catch ensures airtight closure onto a silicone rubber seal.

Wire shelves and a stainless steel desiccant tray are included.



Model	DC1	DC2	DC3	DC4	DC5
Int. Dim. (mm) H x W x D	350 x 300 x 300	460 x 380 x 380	600 x 380 x 380	600 x 500 x 500	900 x 610 x 510
Ext. Dim. (mm) H x W x D	380 x 320 x 320	500 x 400 x 400	620 x 400 x 400	620 x 520 x 520	920 x 650 x 550
Capacity (L)	32	66	90	150	270
Shelves Supplied	2	2	2	4	4
Material (door)	Acrylic				
Material (exterior)	Stainless Steel				
Options	Gas inlet/outlet valves, additional shelves, glass door.				

Serving the Scientific, Medical and Research industries since 1945

LABEC



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au