

Autoclave - Bench Model

Shinva 45 Litre



Serving the Scientific, Medical and Research industries since 1945

LABEC

LABEC offers the **Shinva 45 Bench Top Autoclave**. It is the largest B Class Bench Top autoclave available in Australia.

TGA Approved, for the sterilisation of wrapped and unwrapped solid or hollow instruments. The unit has 4 sterilisation programs, Crystal Display Screen, Auto protection Device and Door Safety Interlock.

The removable racking system is made up of sliding baskets which can be all removed. A vertical racking system is also available as an option.

NATA certification is available and can be supplied at a small cost.

Features include:

- Sterilisation of wrapped & unwrapped solid or hollow items.
- Unit has 4 sterilisation programs.
- Crystal Display Screen.
- Auto Protection Device.
- Door Safety Interlock.



SPECIFICATIONS

Model	Shinva-45B	Shinva-45N
Class of Autoclave	B Class	N Class
Sterilising Function	Solid and Hallow Instruments	Liquids
Internal Dimensions	Ø316 x d618mm	
External Dimensions	w640 x d830 x h550mm	
Capacity (L)	45	
Power (W)	3600	
Power Source	240V / 15AMP	
Controller	LCD Display with adjustable pre-set programs	
Printer	Included	
Weight	120kg	



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au