

# Cold Plate

(+5°C to -20°C)

**LABEC Cold Plate** is used for insect paralysis testing and wax solidification.

Featuring a large plate area and stepped platform for easy use at bench height. The temperature is settable between +5°C and -20°C.

This cold plate is the ideal choice suitable for all routine laboratory applications.



## SPECIFICATIONS

Model	R-TEC-BLIV
Temperature Range (°C)	-20°C to +5°C
Plate Area (mm)	315 x 380
Overall Dimension (mm)	w710 x d350 x h388
Weight (kg)	30kgs
Power (AMPS   V)	10AMPS   240V   300W

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: [sales@labec.com.au](mailto:sales@labec.com.au)

Ph: 02 9560 2811 • Fax: 02 9560 6131

[www.labec.com.au](http://www.labec.com.au)