



## TISSUE FLOATING BATH&SLIDE WARMER ; MODEL NO.: RD-2018



### ABOUT THE PRODUCT

- 1.The baking sheet part is colored with special technology, which has the advantages of wear resistance and no paint peeling.
- 2.Multi-layer bevel design can process 60 tissue slices at the same time.
- 3.This device uses a medical-grade high-definition color touch screen, all information is clear at a glance and easy to operate
- 4.Precise temperature control by computer, fast and stable temperature rise.
- 5.It takes up to 10 minutes to reach 75 degrees.
- 6.The oven dryer and the slide dryer can be set to be heated separately.
- 7.Timing switch function Set time can be accurate to minutes.



### DETAIL PRESENTATION

# RD-2018

TISSUE FLOATING BATH&SLIDE WARMER  
Focus On Pathpology Since 2012.




Product Name	Tissue floating bath & slide dryer
Model No	RD-2018
Floating Bath	0°C ~ 75°C or 0°C ~ 99°C optional
Dryer	0°C ~ 75°C or 0°C ~ 99°C optional
Oven	0°C ~ 75°C or 0°C ~ 99°C optional
Product size	570mm*430mm*180mm
Water bath size	210mm*170mm*55m
Oven size	240mm*152mm*110mm; 60 tissue slices at the same time
Floating Bath	2L
Voltage supply	AC220V±10%, 50Hz±1Hz
Power	AC220V±10%, 50Hz±1Hz


SWITCH BUTTON



TOUCH SCREEN



HEATING PLATE



GLASS TRAY

