

TISSUE FLOATING BATH ; MODEL NO.: XH-1001

ABOUT THE PRODUCT

- 1. microprocessor control and modular, high precision, good stability, easy to operate.
- 2. steel and with white printing coating, good thermal conductivity.
- 3. adjustable temperature, doesn't break the tissue.
- 4. the transparent glass lid, keeps heat evenly, reduces heat loss.

DETAIL PRESENTATION

XH-1001

TISSUE FLOATING BATH



Product Name	Tissue Floating Bath
Outside Dimensions	Lab/Hospital/Histology/Pathology
Product Size(mm)	Diameter 270mm, H 210mm
Bath Dimensions	Diameter 210mm, Depth 65mm
Bath Temperature Range	Ambient-75°C
Capacity	4L
Voltage	110V, 60HZ / 220V, 50HZ
Power	400W
Relative Humidity	less Than 85%
Net Weight	5LBS
Dimensions(package)	330*330*175mm
Gross Weight	8LBS
Customized Support	OEM, ODM
Certificate	CE

