

Soxhlet Water Bath

(Ambient +5°C to 100°C)

LABEC Soxhlet Water Bath

is manufactured with high quality materials with 304 grade stainless steel interior for easy maintenance and excellent durability.

Using an internal element connected to a digital controller the bath temperature is maintained. At a constant and accurate temperature. A manual reset switch is fitted to turn off the elements should the water boil dry.

The baths are supplied with concentric rings which can be adjusted to suit different flask sizes and is fitted with a magnetic drain pipe.



SPECIFICATIONS

Model	F-SX6D
Internal Dimension (mm) H x W x D	120 x 945 x 168
External Dimension (mm) H x W x D	215 x 1100 x 250
Type	6 hole opening type
Temp. Min	Ambient +5°C
Temp. Max	100°C
Temperature Accuracy	+/- 0.5°C
Temperature Controller	P.I.D controller
Thermocouple	Pt 100 ohms
Heater (sheath)	1500 W
Capacity (L)	19 litres
Material Interior	Stainless steel plate
Material Exterior	Steel plate with powder coated finish
Weight (kg)	15