

# Ultrasonic Bath

Small and large capacity

**LABEC** FXP range of **Ultrasonic Cleaners** is a single chamber device used for Medical, Dental, Laboratory purposes where all forms of contamination are lightly deposited.

Features include a single tank system with rounded corner reflex design for easy cleaning and the polyester coated transducers to prevent moisture contamination and maintain high efficiency levels. All units incorporate Lateral Oscillation Technology for uniform cleaning results. The internal basket system allows the lid to seal the working solution avoiding evaporation and reducing ultrasonic acoustic levels. The chamber is pressed 304 stainless steel and 1mm in thickness.

**NOTE:** Lid and basket included.

Sonic energy is provided to the chamber by piezoelectric transducers bonded to the tank bottom. Maximum operating temperature should not exceed 60°C so as to maintain reliability and maximize efficiency.

The clear and easy design with compact and quiet operation makes these products the market leaders.

**OPTIONS:**

- Digital Electronic Timer.
- Heated Units.

Small Capacity



Large Capacity

Serving the Scientific, Medical and Research industries since 1945

**LABEC**

**LABEC**

Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd  
email: sales@labec.com.au  
Ph: 02 9560 2811 • Fax: 02 9560 6131  
www.labec.com.au

# SPECIFICATIONS

## Ultrasonic Bath Small and large capacity

Model	Capacity	Timer	Heating	External Dim (mm)	Internal Dims (mm)	Basket Dims (mm)
FXP8M	1.7L	Mechanical (30min)	Nil	L225*W165*D240	L150*W140*D100	L120*W100*D50
FXP8D		Digital (120min)				
FXP10M	2.7L	Mechanical (30min)	Nil	L315*W165*D240	L240*W140*D100	L210*W100*D50
FXP10D		Digital (120min)				
FXP10MH		Mechanical (30min)	100W			
FXP10DH		Digital (120min)				
FXP12M	5.6L	Mechanical (30min)	Nil	L375*W180*D270	L295*W150*H150	L250*W110*D110
FXP12D		Digital (120min)				
FXP12MH		Mechanical (30min)	175W			
FXP12DH		Digital (120min)				
FXP14M	10L	Mechanical (30min)	Nil	L375*W270*D290	L295*W240*H150	L250*W200*D100
FXP14D		Digital (120min)				
FXP14MH		Mechanical (30min)	275W			
FXP14DH		Digital (120min)				
FXP16M	14L	Mechanical (30min)	Nil	L395*W320*D290	L320*W295*D150	L280*W260*H100
FXP16D		Digital (120min)				
FXP16MH		Mechanical (30min)	300W			
FXP16DH		Digital (120min)				
FXP20M	21.5L	Mechanical (30min)	Nil	L575*W325*D290	L495*W295*D150	L450*W250*H100
FXP20D		Digital (120min)				
FXP20MH		Mechanical (30min)	400W			
FXP20DH		Digital (120min)				
FXP40M	61L	Mechanical (30min)	Nil	L690*W535*D370	L600*W495*D200	L550*W450*D150
FXP40D		Digital (120min)				
FXP40MH		Mechanical (30min)	800W			
FXP40DH		Digital (120min)				

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au