

# Ultra Low Temperature Chest Freezer

*Economy (-40°C/-60°C/-86°C/-135°C/-150°C)*

**LABEC** offer a range of **Economy Ultra Low Temperature Chest Freezers** ideal for the storage of non-critical product. Suitable for industrial and scientific research applications, and general purpose use where low temperatures are required.

**Features include:**

- Digital Control system with over and under temperature alarms.
- Anti corrosion stainless steel liner.
- Stable temperature with double-sealed lid gasket.
- CFC-free refrigerant system.
- Washable condenser filter.
- Powder coat steel finish.



Model H-DW-40W118



Model H-DW-40W58



Model H-DW-40W358

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd  
email: [sales@labec.com.au](mailto:sales@labec.com.au)  
Ph: 02 9560 2811 • Fax: 02 9560 6131  
[www.labec.com.au](http://www.labec.com.au)

## SPECIFICATIONS

# Ultra Low Temperature Chest Freezer

*Economy (-40°C/-60°C/-86°C/-135°C/-150°C)*

Model	External Dims	Internal Dims	Temperature	Capacity
	W*D*H (MM)	W*D*H (MM)	(°C)	(L)
H-DW-40W58	670*765*1025	410*410*385	-15 TO -40	58
H-DW-40W118	1270*795*990	540*400*560	-15 TO -40	118
H-DW-40W258	1510*895*1090	780*500*660	-15 TO -40	258
H-DW-40W358	1810*895*1090	1080*500*660	-15 TO -40	358
H-DW-40W458	2115*895*1095	1380*500*660	-15 TO -40	458
H-DW-60W58	670*765*1025	410*410*385	-25 TO -60	58
H-DW-60W118	1270*795*990	540*400*560	-25 TO -60	118
H-DW-60W258	1510*895*1090	780*500*660	-25 TO -60	258
H-DW-60W358	1810*895*1090	1080*500*660	-25 TO -60	358
H-DW-60W458	2115*895*1095	1380*500*660	-25 TO -60	458
H-DW-86W58	670*765*1025	410*410*385	-40 TO -86	58
H-DW-86W118	1270*795*990	540*400*560	-40 TO -86	118
H-DW-86W258	1510*895*1090	780*500*660	-40 TO -86	258
H-DW-86W358	1810*895*1090	1080*500*660	-40 TO -86	358
H-DW-86W458	2115*895*1095	1380*500*660	-40 TO -86	458
H-DW-135W58	670*765*1025	410*410*385	-100 TO -135	58
H-DW-135W118	1270*795*990	540*400*560	-100 TO -135	118
H-DW-135W258	1510*895*1090	780*500*660	-100 TO -135	258
H-DW-150W258	1510*895*1090	780*500*660	-100 TO -150	258

Serving the Scientific, Medical and Research industries since 1945

**LABEC**



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd

email: sales@labec.com.au

Ph: 02 9560 2811 • Fax: 02 9560 6131

www.labec.com.au