

Sand Bath (up to 350°C)



Serving the Scientific, Medical and Research industries since 1945

LABEC Electric Sand Baths are available in both analogue and digital control. These are made with robust construction for continuous performance, with an aluminium heating plate topped by a stainless steel sand surround. These are suitable for all routine laboratory applications. They are easy to use and durable, and are available in three sizes.



Digital Model LA-SH-9BSB

SPECIFICATIONS

Model	LA-SH-5SB	LA-SH-8SB	LA-SH-9SB	LA-SH-5BSB	LA-SH-8BSB	LA-SH-9BSB
Max Temp.	320°C			350°C		
Shell Material	Painted Mild Steel					
Sand Surround	Aluminum tensile					
Plate Material	Stainless Steel surround, 50mm deep					
Heating power (W)	1200	1800	3000	1200	1800	3000
Temperature Controller	Analogue Dial			Digital		
Temperature Accuracy	±10%			±1%		
Plate Temp. Uniformity	±15%			±10%		
External Dims (mm)	W340 x D370 x H270	W420 x D500 x H270	W470 x D650 x H270	W340 x D370 x H270	W420 x D500 x H270	W470 x D650 x H270
Sand Surround (mm)	W300 x D300 x H50	W350 x D450 x H50	W400 x D600 x H50	W300 x D300 x H50	W350 x D450 x H50	W400 x D600 x H50
Max Plate Load (kg)	20	30	40	20	30	40
Power Supply	240V / 10Amp		240V / 15Amp	240V / 10Amp		240V / 15Amp
Weight (kg)	7	10	15	7	10	15



Laboratory Equipment Pty Ltd

Laboratory Equipment Pty Ltd
 email: sales@labec.com.au
 Ph: 02 9560 2811 • Fax: 02 9560 6131
 www.labec.com.au