



Y-LS-400, Floor Sieve Shaker

- Standard sieve-shaker model, with both reciprocating and orbital shaking at same time
- With a continuously striking device on the top for helping sieving jobs well.
- No sieves included
- Please confirm Tapping should be 156 tpm or 78 tpm
- Available to use the sieves in diameter of 20cm.

| Model | Y-LS-400 |
|---------------------|-----------------------------|
| Speed (EPM) | 246 (50Hz) / 291 (60Hz) epm |
| TaPPing | 156 / 78 tpm |
| Time Setting | ON/ 0 ~ 60min |
| Dimensions | 400x750x950mm |
| Sieve Size | φ200x60mm |
| Weight | 185kg |
| Power | AC110/220V, 50/60Hz |