


TISSUE FLOATING BATH ; MODEL NO.: RD-80



ABOUT THE PRODUCT

- 1.LED touching screen allows user to see all temperatures without scrolling, intuitive control panel.
- 2.Fully programmable ON/OFF control allows automatic system start and stop anytime.
- 3.Low-voltage, safe and bright LED lamp, and allows for comfortable processing of specimens.
- 4.Large well drained mould oven,8 holes can easy to open drainage trays under all work areas.



DETAIL PRESENTATION

RD-80

TISSUE FLOATING BATH

Focus On Pathology Since 2012.



Model:	RD-80
Product size:	327mm*340mm*166mm(L*W*H)
Product weight:	8KGS
Temp range:	0°C ~ 70°C, precise 1°C;
Dimensions of bath bowl:	210mm*168mm*55mm(L*W*H)
Capacity:	2.5L
Voltage and frequency:	AC220V±10%, 50Hz±1Hz; 110V accepted custom made
Light at the bottom:	Yes
Package size:	400mm*390mm*240mm(L*W*H); 9.5KGS
Certification:	CE
Warranty:	1-Year
Custom-made:	Accept

SWITCH BUTTON



TOUCH SCREEN



GLASS TRAY



HEAT SINKS

