

SI-400, Vortex Mixer



SI-100N

Features:

- Compact, robust, and durable design for efficiency.
- Variable speed control.
- Continuous operation or automatically press start.
- Fast thorough, and quiet vortex mixing action, but vibration free.
- Using high quality motor to achieve the performance.
- Option 1: with a platform head suitable for mixing the contents of several test tubes simultaneously or a flask.
Option 2: with a single cup head for mixing the content of one tube.
- Quality certificate ISO 9001, ISO 13485.

Model	SI-100N
Max. Speed	3000rpm
Orbit	5.0mm
Top Head Dimen.	Rubber platform head with Φ 80mm Rubber cup head with Φ 26mm
Dimensions(WxDxH)	123x163x170mm
Weight net/gross	3.7/4.2kg