

## Table Top Low Speed Centrifuge

### Introduction:

BKC-TL4CII centrifuge is a table top low speed centrifuge with low speed. It is suitable for the separation of blood samples in clinical laboratory, laboratory, beauty salon, dental hospital, etc.

### Features:

- \* With dedicated CGF program.
- \* Small size, safe and reliable, easy to use.
- \* The body is made of high quality plastic material, explosion-proof protection.
- \* 10 kinds of lifting speed to choose from, the lifting speed is fast.
- \* Brushless DC motor, stable operation, small vibration, low noise.
- \* Push-button programming design, real-time display of operating parameters, more intuitive operation.
- \* Can set the centrifugal speed, time, speed up and down gear, and check the centrifugal force at any time.
- \* Ultra-low noise, microcomputer control, intelligent interface, fault display, door cover, overspeed and other automatic protection functions.



### Technical Parameters:

Model	BKC-TL4CII
Max.Speed	4000rpm
Max.RCF	1880×g
Max.Capacity	8*20ml
Speed Precision	±20rpm
Time Range	1s~999min
Power Supply	AC110V/220V 50/60Hz
Consumption	65W
Noise	≤55dB(A)
External Size(L*W*H)	370*270*210mm
Package Size(L*W*H)	370*310*250mm
Net/Gross Weight	7.5kg/8.5kg

### Rotor:

No.	Rotor	Capacity	Max. Speed(rpm)	Max.RCF(×g)
No.1	Angle Rotor	8*20ml	4000	1880
No.2		8*15ml (Suitable for 5/7/10ml vacuum tube)	4000	1840
No.3		12*15ml (Suitable for 5/7/10ml vacuum tube)	4000	1790

### Program Selection:

Program	CGF				A-PRF	I-PRF	PGRF
Speed(rpm)	2600	2300	2600	2850	2450	2850	3100
Time(min)	2	4	4	3	4	4	3